

# **PRATHAM<sup>®</sup> SYPHILIS / TP ELISA**

*ELISA for the detection of Antibodies to Syphilis in human serum / plasma.*

## **INTENDED USE**

**PRATHAM<sup>®</sup> Syphilis / TP ELISA** Test is intended to be used for the detection of all subtypes (IgG, IgA, IgM) of antibodies to Syphilis in human serum or plasma.

## **INTRODUCTION**

*Treponema pallidum* (TP) is the causative agent of venereal disease Syphilis. TP is a spirochete bacterium with an outer envelope and a cytoplasmic membrane. After infection, host forms non –treponemalantipodial antibodies to the lipodial material released from the damaged host cells as well as treponema specific antibodies. Diagnosis of syphilis depends on the correlation of clinical data with the non-treponemal and treponemalassays. Non-treponemal tests (VDRL, RPR, etc) are generally used for screening and treponemal tests (TPHA, FTA-ABS) are used as confirmatory tests. Treponema antibody ELISA tests are gaining importance as screening and conformity tests, as they detect the presence of antibodies specific to *Treponema pallidum*.

## **TEST PRINCIPAL**

**PRATHAM<sup>®</sup> Syphilis / TP ELISA** micro well strips are coated with recombinant antigens of syphilis. Samples along with positive and negative controls are added in the coated wells and incubated. The wells are washed to remove unbound components and syphilis recombinant antigens conjugated to Horse Reddish Peroxides are added. After a short incubation the wells are washed again and bound enzyme is detected by adding substrate. The reaction is stopped after specified time with acid and absorbance is determined for each well at 450nm with an ELISA reader. The cut-off value is calculated by the given formula and absorbance of all the wells are compared with the cut-off value. Any sample having absorbance more than the cut-off value is considered reactive.

## **KIT COMPONENTS**

Component	Description of Reagent	Presentation
Coated Microwells	96 Microwell (12x8) are coated with recombinant Syphilis Antigens.	1
Positive Control	Inactivated and stabilized human serum reactive for Syphilis with preservatives.	1x1.5 ml
Negative Control	Inactivated and stabilized human serum non-reactive for HIV-1, Syphilis HBsAG and HCV.	1x1.5 ml
Sample dilution buffer	Buffered solution containing stabilizing proteins and preservatives.	1x6 ml
Conjugate	Recombinant Syphilis-HRP conjugate (21X). To be diluted 21 times with conjugate dilution buffer.	1x0.8 ml
Conjugate dilution buffer	Buffered solution containing proteins and preservative.	1x12 ml
Substrate	Solution containing Tetramethibenzidine (TMB) and hydrogen peroxide. Ready to use.	1x10 ml
Wash buffer	Buffer containing surfactants (40X). To be diluted 40 times with distilled or deionized water.	1x25 ml
Stop solution	Diluted sulphuric acid.	1x10 ml
Plate sealer	Adhesive paper sheet for covering plate on incubation time	2 Nos.
Pack insert	Instruction for use	1 Nos.

## **MATERIAL REQUIRED BUT NOT PROVIDED**

Microplate reader capable of measuring absorbance at 450nm, Pipettes and Pipettes tips, Deionized or distilled water, 1X PBS, Automated microplate washer, Syphilis Reader, Disinfectant, Timer, Gloves, Biohazard Waste Container etc.

## STORAGE AND STABILITY

**PRATHAM® Syphilis / TP ELISA** test components are stable up to expiry date indicated on the component label/box label. Syphilis ELISA kit needs to be stored at 2-8°C.

## SPECIMEN COLLECTION AND STORAGE

1. No prior preparation of the patient is required.
2. Collect blood specimen by venipuncture according to the standard procedure.
3. Specimen should be free of particulate matter and microbial contamination.
4. Preferably use fresh sample. However, specimen can be stored refrigerated for 24 hours. For long storage, freeze at - 20°C or below. Specimen should not be frozen and thawed repeatedly. Maximum of two freeze/thaw cycles are allowed.
5. Do not use heat inactivated specimen.
6. Specimen containing precipitate or particulate matter should be clarified by centrifugation prior to use.
7. No not use turbid, lipaemic, haemolysed, clotted or contaminated specimen.

## WARNING AND PRECAUTIONS

- Do not use expired reagents.
- Do not mix reagents from different lots within a given test run.
- Before use, allow reagents to reach room temperature (+18-30°C).
- Carefully reconstitute or dilute the reagents avoiding any contamination.
- Do not carry out the test in the presence of reactive vapors (acid, alkaline, aldehyde vapors) or dust that could alter the enzyme activity of the conjugate.
- Use glassware thoroughly washed and rinsed with deionized water or preferably, disposable material.
- Do not allow the micro plate to dry between the end of the washings operation and the reagent distribution.
- The enzyme reaction is very sensitive to metal ions. Consequently, do not allow any metal element to come into contact with the various conjugate or substrate solutions.
- Use a new pipette tip for each sample.
- Washing the micro plate is a critical step in the procedure: Follow the recommended number of washings cycles and make sure that all wells are completely filled and then completely emptied. Incorrect washings may lead to inaccurate results.
- Inadequate removal of residual wash buffer can cause inconsistent color development. Microwell strips should be tapped and blotted on absorbent paper or towels to minimize residual wash buffer.
- Never use the same container to distribute conjugate and development solution.
- Check the pipettes and other equipment for accuracy and correct operations.
- Do not change the essay procedure.

### Preparation of Reagents prior to use:

**Wash Buffer Preparation:** Dilute wash buffer 40 times for example add 5 ml concentrated buffer to 195 ml distilled or deionized water.

**Conjugate Preparation:** Dilute conjugate 21 times as shown below:

Reagent (Strip)	1	2	3	4	5	6	7	8	9	10	11	12
Conjugate (µl)	50	100	150	200	250	300	350	400	450	500	550	600
Conjugate dilution buffer (ml)	1	2	3	4	5	6	7	8	9	10	11	12

## TEST PROCEDURE

1. Bring all the reagents and specimen to room temperature before use.
2. Take out required number of strips and immediately close the pouch.
3. Prepare data sheet indicating the location of controls and specimen.
4. Use controls in duplicate.
5. Leave well A1 as substrate control.
6. Add 50µl sample dilution buffer in each well except A1.
7. Add 100µl sample or controls in separate well except well A1.
8. Apply plate sealer and incubate for 60 minutes at 37°C.

9. Wash each well by filling approximately 350µl diluted wash buffer and aspirating/flicking off six times Blot dry.
10. Add 100µl diluted conjugate in each well except A1 and incubate for 30 minutes at 37°C.
11. Wash six times as in step 9 and blot dry.
12. Add 100µl substrate in each well including A1 and incubate at room temperature away from light for 15 minutes.
13. Stop reaction by adding 100µl stop solution. The stop solution should be added in the same sequence as substrate addition.
14. Blank ELISA reader with well "A-1".
15. Read the absorbance at 450nm with 630nm or above as reference within 30 minutes of stopping the reaction.

### Quality Control

**Negative Control:** The individual absorbance value of negative controls should be less than 0.1.

**Positive Control:** The individual absorbance value of positive controls should be more than 1.0.

### Calculation of Cut of Vluue

The cut off value =  $0.2 + NC(Avg)$

### Set up the cut-off value

#### Example

NC	Absorbance
A1	0.022
B1	0.026

**Avg. of NC=**  $(0.022+0.026)/2= 0.024$

cov =  $0.2+0.024=0.224$

## INTERPRETATION OF RESULTS

1. Samples with absorbance value equal to or less than the cut-off value is considered non-reactive by Syphilis ELISA kit and are considered negative for syphilis antibodies.
2. Samples with absorbance value greater than cut-off value are considered reactive by Syphilis ELISA kit. The reactive samples should be retested in duplicate. Initially reactive sample that do not react in either of duplicate are considered negative for antibodies to syphilis. Initially reactive sample that reacts in either or both duplicates are considered repeatedly reactive.
3. If a sample is repeatedly reactive the probability of antibodies to syphilis are high, especially with patients at high risk or high absorbance values. Such samples should be retested with supplemental tests.

## INTERNAL QUALITY CONTROL

1. Positive and negative controls should be included in each test batch.
2. Acceptable performance is indicated when a uniform milky suspension with no agglutination is observed with the Dengue NS1 Negative Control and agglutination with large aggregates is observed with the Dengue NS1 Positive Control.

## PERFORMANCE CHARACTERISTICS

The kit has been evaluated with the 50 known panel of Dengue NS1 Ag positive and 100 negative samples. Following is the in-house evaluation.

Description	Negative	Positive	Total
Syphilis known Negative Specimen	100	0	100
Syphilis known Positive Specimen	0	50	50

Assay Sensitivity : 100%

Assay Specificity : 100%

## LIMITATION OF THE TEST

1. Though Syphilis ELISA is a reliable screening assay, it should not be used as a sole criterion for diagnosis of syphilis infection. Reactive sample should be retested with confirmatory assays.
2. Absence of syphilis antibody does not indicate that an individual is absolutely free of syphilis infection.
3. Since various tests for syphilis differ in their performance characteristics and antigen composition, their reactivity patterns may differ.
4. Testing of pooled samples is not recommended.
5. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

## DISPOSAL

Consider all test run with human specimen as potentially infectious and discard using standard biosafety practices.

## DISCLAIMER:











Whilst every precaution has been taken to ensure the diagnostic ability and accuracy of this product the product is used outside of the control of the Manufacturer and Distributor and the result may accordingly be affected by environmental factors and / or user error. A person who is the subject of the diagnosis should consult a doctor for further confirmation of the result.

The manufacturer and distributors of this product shall not be liable for any losses, liability, claims, costs or damages whether direct or indirect or consequential arising out of or related to an incorrect diagnosis, whether positive or negative in the use of this product.

## REFERENCES

1. Aral R. Marx. Crack, sex and STD, Sexually Transmitted Diseases, 1991;18:92-101.
2. J.N. Wasserheit. Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases, Sexually Transmitted Diseases 1992;19:61-77.
3. Johnson Phillip C. Testing for Syphilis, Dermatological Clinic 1994;12 Jan:9-17.
4. Clarie M. Fraser. Complete genome sequence of Treponema Pallidum the Syphilis spirochete, Science 1998;281 July: 375-381.
5. Center for Disease Control. Recommendations for diagnosing and treating Syphilis in HIV-infected patients, MMWR Morb. Mortal Wkly Rep. 1998; 37:607.

## SYMBOLS

 Read instructions for use	 Name of Manufacturer	 For single use only
 No. of test	 Expiry Date of Kit.	 Date of manufacturing of IVD Kit
 In-vitro diagnostic use	 Keep away from Sunlight	 Reference Catalogue Number
 Storage Condition		